

Calibration Kit

For NO and NO₂



When the smallest thing matters

Specifications

Everything you need for NO and NO₂ gas calibration

The SLE NO and NO₂ Calibration kit, with the economy of disposable gas canisters

- Accurate, reliable and easy gas calibration
- 2 disposable gas canisters each containing 58 litres of calibration gas
- More economical than re-usable cylinders
- A re-usable, interchangeable safety pressure regulator comes as standard in the kit
- A robust and lightweight carrying-case gives a total kit weight of only 3.5 kg offering greater convenience and portability
- No cylinder refill problems or rental charges

Ordering Details

Z4123

Complete kit

Consisting of:

N4123/01

Disposable canister of NO/Nitrogen gas mixture containing a nominal 58 litres of 25 ppm NO. (Certified NO concentration between 20-30 ppm in pure nitrogen)

N4123/02

Disposable canister of NO₂/Air gas mixture containing a nominal 58 litres of 10 ppm NO₂. (Certified NO₂ concentration between 5-15 ppm in pure air)

N4123/03

Reusable / interchangeable gas regulator

N4123/04

Lightweight carrying case - 19 (W) x 56 (H) x 10 (D) cm



SLE Limited.

Twin Bridges Business Park, 232 Selsdon Road,
South Croydon Surrey CR2 6PL UK

Telephone: +44 (0)20 8681 1414 • Fax: +44 (0)20 8649 8570

E-mail: sales@sle.co.uk • Web: www.sle.co.uk